

Dräger DrugCheck[®] 3000

Oral Fluid POCT Drug Detection Device

- Quick & Easy to Use
- Compact, Portable Device
- No Power Required
- Extremely Hygienic
- Proven Dräger Quality
- Measures THC at 15 ng/ml



The Dräger DrugCheck 3000 is a compact oral fluid POCT drug detection device that delivers results within minutes.

Benefits

- Hygienic use and disposal
- Compact, pocket-size device
- Easy to read results
- Virtually impossible to adulterate
- Use anywhere, even remote areas without power
- Sample collection and results reading, all in one device
- Results are ready in 5 - 7 minutes
- Adjustable cut-off for THC: 15 or 25 ng/mL (fast mode)

Usage

The Dräger DrugCheck 3000 consists of two components: an absorbent collector for obtaining an oral fluid sample, and a test cassette for the analysis. The test cassette contains the buffer liquid and a window with two test strips, which displays the control and test lines. The drug test itself is performed in three easy steps: swab the saliva sample, shake the test kit, wait for a brief incubation period, and then start the test. A colour indicator on the swab disappears as soon as it has absorbed enough oral fluid for a test.


As soon as the control lines appear in the window, you can read the results. If the test result is negative, a line will appear alongside the respective drug class. This means that none of the target drug was detected in the sample. If a line fails to appear next to one of the drug classes, then the result for this drug is non-negative.

Drugs & Pharmaceuticals Detected

AMP	Amphetamine	35 ng/ml
MET	Methamphetamine	35 ng/ml
COC	Cocaine/Benzoyllecgonine	20 ng/ml
OPI	Morphine	20 ng/ml
THC	Δ^9 -tetrahydrocannabinol	15 or 25 ng/ml



Components



Specifications

Packaging:	Box of 20 individually sealed tests
Dimensions:	32 x 111 x 57 mm (1.26 x 4.37 x 2.2 in)
Weight:	30 grams
Operating Temp:	5°C to 30°C (41°F to 86°F)
Measurement Time:	< 5 min (fast mode) < 7 min (standard mode)

SBL Testing Technologies Inc. (SBL) offers the leading oral fluid drug testing and breath alcohol technologies available globally. SBL also provides an unparalleled platform of comprehensive support services to ensure our client's success with these technologies. From legal and policy to training and technical, SBL's experienced team of subject matter experts will ensure success.

At SBL, we understand the sensitive nature of workplace drug & alcohol testing. Our experience and expertise will ensure that your company is ready to integrate and optimize all the advantages of oral fluid drug testing. With a corporate philosophy of unwavering commitment to service, SBL is your partner to enhance your drug & alcohol program: seamlessly and cost-effectively while successfully addressing the dynamic challenges faced by corporate drug & alcohol programs today in Canada.

In Canada, the DrugCheck 3000 is for **FORENSIC USE ONLY**

403 966 5450
info@sbltt.com
www.sbltt.com

